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**IDENTIFICATION OF FOOT-AND-MOUTH DISEASE  
VIA ANTIBODIES**

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## CONTENTS

	Page
Identification of foot-and-mouth disease VIA antibodies in animal sera: production of reagents, performance of tests and interpretation of results. .	19
— Introduction . . . . .	21
— Materials and Methods. . . . .	22
Preparation of the VIA antigen . . . . .	22
Inactivation of the suspension containing the VIA antigen. . . . .	23
Control serum . . . . .	23
Preparation of the 2% agar . . . . .	24
Glycine buffer. . . . .	24
Preparation of plates . . . . .	24
Titration of the reagents . . . . .	24
Determination of anti-VIA antibodies in animal sera . . . . .	25
Quantification of VIA antibodies by IDAG titration. . . . .	25
— Results . . . . .	26
Titration of the VIA antigen and the control serum . . . . .	26
Analysis of sera to identify VIA antibodies. . . . .	26
Titration of sera . . . . .	26
— Remarks . . . . .	30
— References. . . . .	31

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---

IDENTIFICATION OF FOOT-AND-MOUTH DISEASE  
VIA ANTIBODIES IN ANIMAL SERA:  
PRODUCTION OF REAGENTS, PERFORMANCE OF  
TESTS AND INTERPRETATION OF RESULTS

---

## INTRODUCTION

A virus-infection-associated (VIA) antigen was identified in foot-and-mouth disease (FMD) virus suspensions (4). The antigen appeared to be an inactive form of the ribonucleic acid (RNA) viral polymerase (8). Antibodies to the VIA antigen were detected by the immune diffusion agar gel (IDAG) test of sera from cattle that had been affected by FMD (6) or had been inoculated with either attenuated live viruses or formalin-inactivated vaccines (1). However, sera from animals first vaccinated with vaccines inactivated with inactivants of a first-order kinetics showed no anti-VIA antibodies detectable by IDAG (1, 6). However, when these animals were inoculated with FMD virus, many of them developed VIA antibodies without showing signs of the disease (1). Recent studies (5, 7) have demonstrated that repeated application of inactivated vaccines induce in cattle the appearance of anti-VIA antibodies, albeit of low intensity and short duration.

The VIA antibodies and antigen are characteristically FMD specific, but not virus type specific (4). Identification of these antibodies therefore serves to detect infection, but does not specify the type of active virus.

From the foregoing, it is deduced that identification of anti-VIA antibodies by means of IDAG tests provides a valuable aid in the study of the prevalence of the FMD virus in the animal populations susceptible to the disease and in the selection of FMD-free herds in endemic areas, which herds could provide animals for use in immunological testing or for exportation. Correct interpretation of the results obtained in the IDAG tests, in conjunction with an adequate sampling of the population to be examined, enables to ascertain whether a herd

has or has not had recent contact with the FMD virus. In case of doubt, tests can be undertaken to isolate virus from the oesophageal-pharyngeal (OP) fluid (9).

This technical report describes the methodology utilized by the Pan American Foot-and-Mouth Disease Center (PAFMDC) for the production of VIA antigen and for the quantification of the VIA reagents. The procedure followed for sera analysis by means of IDAG and the interpretation of the results are also explained.

## MATERIALS AND METHODS

### *Preparation of the VIA antigen*

BHK<sub>21</sub>, Clon 13 cell suspensions, cultivated in roller bottles for 48-72 hours with 100 ml of modified Eagle medium not supplemented with serum, are inoculated with 10<sup>5</sup> infective dose 50% (ID<sub>50</sub>) of virus and kept at 37°C. When the monolayer shows a marked cytopathic effect and the cells begin to detach—this occurs approximately 18 hours post inoculation—5 ml of chloroform are added to each bottle. Then the bottles are again shaken, cooled to -20°C, thawed out, and clarified by centrifugation at 1000 g for 10 minutes.

0.05% (w/v) of DEAE-Sephadex A<sub>50</sub> (0.5 g per liter of suspension) is then added to the clarified virulent suspensions. The suspension is then stirred lightly at 4°C for 48 hours, after which the Sephadex is collected by straining the mixture through a 50 ml syringe having a glass wool in the bottom kept in place by a fine metal screen.

Once the Sephadex column is prepared, it is washed with buffer A (0.15M NaCl and 0.02M Tris at pH 7.6) until the phenol red is eluted from the Sephadex. Then a very slow passage of buffer B (1M NaCl and 0.02M Tris at

pH 7.6) is made, at approximately 150 ml per gram of Sephadex. This buffer, which elutes the VIA, is collected in 30 to 50 ml fractions to which is added a saturated (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> solution until turbidity is produced by the precipitation of the proteins.

The precipitate is removed by centrifugation at 8000 g for 20 minutes at 4°C and resuspended in buffer A at the proportion of 5 ml per 1000 ml of virulent virus suspension. Later is clarified by centrifugation at 3000 g for 20 minutes at 4°C. The supernatant contains the VIA antigen.

### *Inactivation of the suspension containing the VIA antigen*

The VIA antigen obtained by means of the aforementioned procedure may contain active virus. It is rendered innocuous through inactivation with 0.05M of binary ethyleneimine (BEI) (2) at 30°C for 12 hours. The innocuity of the inactivated virus suspension is tested by inoculating a roller bottle containing Clon 13 BHK<sub>21</sub> cells with 1 ml of VIA antigen and conducting two consecutive passages at 48 hours post inoculation. The VIA antigen is considered free of active virus when neither passage discloses a cytopathic effect in the cell monolayer, nor is virus detected by the complement fixation test. Then 0.02% of sodium azide is added to the inactivated VIA antigen and it is stored at 4°C.

### *Control serum*

The control serum is a mixture of hyperimmune O, A and C sera obtained from guinea pigs infected and hyper-immunized at 6 weeks post infection with 0.2 ml of live virus administered twice with a one-week interval. Blood samples are taken from the animals one week after the second hyperimmunization and the serum is inactivated for 30 minutes at 56°C. The control serum innocuity test

is conducted in the same manner as indicated for the VIA antigen.

#### *Preparation of the 2% agar*

20 grams of noble agar (Difco) are dissolved in 1100 ml of demineralized water and sterilized for 10 minutes at 15 lbs pressure. The solution is then poured into a tray and left to solidify. It is then cut into approximately 2-cm cubes that are placed in a running waterbath for 24 hours. Later the cubes are washed three times again with demineralized water and kept in the water at 4°C until utilized.

#### *Glycine buffer*

The buffer is composed of 1M glycine, 0.057M sodium diethylbarbiturate, and 1% sodium azide having a pH of 7.8 adjusted with a 1N solution of HCl.

#### *Preparation of plates*

The 2% agar is melted and then mixed with equal parts of the glycine buffer. Then 16 ml of the 1% agar mixture are poured into plastic Petri plates measuring 90 mm in diameter. The plates are placed on a levelled surface and the temperature is kept between 20 and 30°C. After the agar has jelled, wells are made by means of a seven-punch mold. In the center of the mold there is a well punch, around which are symmetrically arranged six punches, each measuring 4-mm in external diameter and 2-mm between them. Four sets of wells are made on each plate. Light suction is used to remove the agar cut from the wells.

#### *Titration of the reagents*

Prior to use each batch of VIA antigen or control serum is titered in IDAG (Figure 1). For this purpose the antigen is diluted 1:1, 1:2, 1:4, 1:8, 1:16 and 1:32 in

glycine buffer and is placed in the six outlying wells of the 4 molds on a Petri plate. The control serum is also diluted 1:1, 1:2, 1:4 and 1:8 in glycine buffer and distributed in the central wells. Each of the cavities should be filled to the brim with the corresponding dilution. A good distribution of the reagents can be obtained by means of Pasteur pipettes and rubber bulbs.

The plates are kept at room temperature on a levelled surface. If the work area is dry, a humidified environment is recommended. The titration reading is taken after 1 to 3 days. The optimum dilution for using each reagent is the highest one that yields a clear reaction.

#### *Determination of anti-VIA antibodies in animal sera*

Petri plates are prepared and used as in the titration of the reagents. 16 sera can then be analyzed on each plate (see Figure 2).

The VIA antigen in the suitable dilution is put in the central wells. The control serum, likewise in the pre-established dilution, is placed in two opposite wells on the edge of each mold. The sera to be tested are distributed in the four remaining wells.

The remaining details of the technique are as indicated in the titration of the reagents.

The reactions are observed between 1 and 3 days. The test will be regarded as valid only when the controls yield clear precipitation bands.

The sera analyzed in Figure 2 show a reaction that is interpreted as indicated in Table 1.

#### *Quantification of VIA antibodies by IDAG titration*

The preparation of the plates, the distribution of the VIA antigen and of the control sera, and the reading, are as described above. The sera that in the previous test have

yielded an intense reaction (similar to that of the control serum) can be titered in IDAG. In this regard, the 1:2, 1:4, 1:8 and 1:16 dilutions are distributed in the four cavities intended for the unknown sera on each mold. Four sera can therefore be titered on a single plate.

The titer of the serum is the highest dilution that yields a positive reaction.

## RESULTS

### *Titration of the VIA antigen and the control serum*

Figure 1 shows the reactions obtained from titration of a VIA antigen and a positive VIA control serum. This test produces a clear reaction of the antigen and the serum up to the 1:16 and 1:4 dilutions, respectively. Nevertheless, the dilutions that are used in routine tests to ensure good reactions are 1:6 for the antigen and 1:2 for the serum.

### *Analysis of sera to identify VIA antibodies*

Figure 2 indicates the results obtained from IDAG analysis of 16 cattle sera. Table 1 shows the interpretation of the reactions of the sera examined in Figure 2. In some instances, some sera yield non-identity reactions; that is, precipitation lines cross the positive serum band (Figure 2, serum No. 4); this indicates that such a reaction is not due to antibodies against VIA. The sera that do not yield identity reactions are regarded as negative. Non-identity reactions have been observed in positive and negative VIA sera.

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Figure 3 shows the scheme for titration of several sera that yielded an intense positive reaction in Figure 2 (sera 5, 6, 12 and 14). Their titers are, respectively, 1:4, 1:8, 1:16 and 1:8.


FIGURE 1. *Titration of VIA antigen and control serum*

- 1, 2, 3, 4, 5 & 6 = Antigen dilutions 1:1, 1:2, 1:4, 1:8, 1:16 and 1:32, respectively.  
a, b, c & d = Serum dilutions 1:1, 1:2, 1:4 and 1:8, respectively.


FIGURE 2. Analysis of cattle sera to identify anti-VIA antibodies by the immune diffusion agar gel test

- a = VIA antigen in 1:6 dilution titered in Figure 1.  
 s = Control serum (VIA positive) in 1:4 dilution titered in Figure 1.  
 1 to 16 = Cattle sera under analysis.

TABLE 1. Interpretation of the reaction provided by the cattle sera examined in Figure 2

Sera No.	Reaction intensity	Sera No.	Reaction intensity
1	+	9	-
2	++	10	+
3	-	11	+
4	-	12	++
5	++	13	++
6	++	14	++
7	-	15	++
8	-	16	+

- = Negative. + = Positive. ++ = Intense Positive.


FIGURE 3. Quantification of anti-VIA antibodies by titration of sera having intense reactions

- a = VIA antigen.  
 s = Control serum (VIA positive).  
 1, 2, 3, 4 = Dilutions 1:2, 1:4, 1:8 and 1:16 of each respective serum.

## REMARKS

The foregoing process utilized to prepare and titer the VIA antigen and the positive VIA control serum has been extensively tested at PAFMDC and in most of the South American countries.

In the IDAG test, the sensitivity and reproducibility of the test is dependent on the use of the reagents at optimal dilutions. Correct titration is therefore of utmost importance.

Quantification of the intensity of the sera's reactions in the routine test, together with an adequate sampling of the population to be examined by dividing it into age groups, is of great utility in attaining the correct characterization of the different epidemiological situations due to the nature of the past experience which the animal has had with the FMD virus.

It is recommended that blood samples be taken from the animals three or more months post vaccination.

The application and interpretation of the anti-VIA antibodies in the field is fully discussed in Scientific Monograph No. 6 issued by the PAFMDC (3).

## REFERENCIAS – REFERENCES

1. ALONSO FERNANDEZ, A., AUGÉ DE MELLO, P., GOMES, I., ROSENBERG, F. El uso del antígeno asociado a la infección viral (VIA) en la detección de ganado expuesto al virus de la fiebre aftosa. (The use of virus-infection-associated antigen (VIA) in the detection of cattle exposed to foot-and-mouth disease virus). *Bol. Centr. Panam. Fiebre Aftosa 17-18*: 17-22, 1975.
2. BAHNEMANN, H.G. Binary ethylenimine as an inactivant for foot-and-mouth disease virus and its application for vaccine production. *Arch. Virol.* 47 (1): 47-56, 1975.
3. CENTRO PANAMERICANO DE FIEBRE AFTOSA. El uso de las pruebas del antígeno asociado con la infección por virus (VIA) de la fiebre aftosa. *Ser. Monog. Cient. Técn.* 6, 1980.
4. COWAN, K.M. & GRAVES, J.H. A third antigenic component associated with foot-and-mouth disease infection. *Virology* 30 (3): 528-540, 1966.
5. DAWÉ, P.S. & PINTO, A.A. Antibody responses to type-specific and virus-infection-associated antigens in cattle vaccinated with inactivated polyvalent foot-and-mouth disease virus in North Malawi. *Br. vet. J.* 134: 504-511, 1978.
6. McVICAR, J. & SUTMÖLLER, P. Foot-and-mouth disease: the agar gel diffusion precipitation test for antibody virus-infection-associated (VIA) antigen as a tool for epizootiologic surveys. *Am. J. Epid.* 92 (4): 273-278, 1970.
7. PINTO, A.A. & GARLAND, A.J.M. Immune response to virus-infection-associated (VIA) antigens in cattle repeatedly vaccinated with foot-and-mouth disease virus inactivated by formalin or acetyleneimine. *J. Hyg. (Camb.)* 82 (1): 41-50, 1979.
8. POLATNICK, J., ARLINGHANS, R., GRAVES, J.H., COWAN, K.M. Inhibition of cell-free foot-and-mouth disease virus ribonucleic acid synthetics by antibody. *Virology* 31: 609-615, 1967.
9. SUTMÖLLER, P. & COTTRAL, G.E. Improved techniques for the detection of foot-and-mouth disease virus in carrier cattle. *Arch. ges. Virusforsch.* 21 (2): 170-177, 1967.